

3M™ Petrifilm™ -Staph Express Disks

Product Information

MEDIUM:	Petrifilm™ Staph Express Disks (6492/6493)				
ISO/BSI CERTIFICATE OF REGISTRATION NUMBER:	3M Food Safety is certified to ISO-9001:2008, FM 14552 Manufactured at Brookings, South Dakota, U.S.A.				
DATE OF EXPIRATION/ LOT NUMBER	Expiry and lot number indicated on each package. Lot number indicated on each plate.				
FORMULATION:	Toluidine Blue-O dye and a DNA indicator				
METHOD OF PREPARATION:	DNA indicator and Toluidine Blue-O dye coated onto films. For use, insert into Petrifilm™ Staph Express Count plates. See product instructions for detailed information.				
EFFICACY CHECK:	The disks are tested for the appropriate functioning of the Toluidine Blue-O DNA indicator system. The following test determines the function: <table><thead><tr><th>Enzyme</th><th>Acceptable Result</th></tr></thead><tbody><tr><td>Micrococcal Nuclease</td><td>Pink color development where enzyme is placed within two (2) minutes</td></tr></tbody></table>	Enzyme	Acceptable Result	Micrococcal Nuclease	Pink color development where enzyme is placed within two (2) minutes
Enzyme	Acceptable Result				
Micrococcal Nuclease	Pink color development where enzyme is placed within two (2) minutes				
PACKAGING:	Pack size: 20 individually packaged disks per foil pouch Film grade: Plastic foil laminate Seal integrity check: Delamination seal integrity test performed				
MEASUREMENT SYSTEMS CALIBRATION AND TRACEABILITY:	Incubator temperature. 3M internal calibration. Minimum calibration once per year for all equipment.				
MEDIA QUALITY STATEMENT:	Quality assurance certificate included in package				
SHELF LIFE:	18 months from date of manufacture				
STORAGE CONDITIONS:	Store at temperatures less than or equal to 8°C				
SIGNED:					



Marta Henrickson
Quality Assurance